

MRI Devices

Instruments for Interventional Magnetic Resonanz Imaging (MRI)



SOMATEX[®]
MEDICAL TECHNOLOGIES GMBH

www.somatex.com

MRI Devices

Optimal visibility and maximum strength

Magnetic resonance imaging (MRI) is a future oriented imaging procedure in modern diagnostics and minimal invasive therapy. Its advantages include better visualization of several organs and extremely high detail recognition – without harmful ionizing radiation. Additionally MRI is adopted for minimally invasive treatments of patients. Instruments for MRI have to meet a very high standard.

SOMATEX® has been active for years in the development of instruments made from non-magnetic material for different applications. In accordance with the high demands and new developments in MRI, our products are continually improved or newly designed.

The SOMATEX® MRI Instruments find the optimal balance of MRI visibility, artefact management, and material strength. Thanks to special design and the use of new high quality materials, SOMATEX® is setting high standards in the MRI sector: Minimum artefacts, optimal visibility, and high strength make the SOMATEX® MRI Instruments safe and reliable.


MRI Devices at a glance

- Tumark® MRI
- MRI Tuloc®
- MRI DUO-System
- MRI Biopsy-Handy
- MRI Coaxial Puncture Needle
- MRI Chiba Needle


Tumark® MRI

High visibility

The **Tumark® MRI** allows for precise tissue marking under MRI control. The Tumark® MRI is adapted in shape and size for the application within the MRI gantry. The placement of the marker can easily be performed inside the MRI. The ergonomic handle with slide button furthermore allows for easy, single-handed deployment. The new 3D marker design enables firm anchorage in the tissue with optimal visibility in all positions. Approved 3D marker design ensures firm anchoring in the tissue with optimal visibility in all positions. For the follow-up, the marker will be clearly visible in Magnetic Resonance, ultrasound, stereotaxy or X-Ray.



Improved visibility for ultrasound



Application

- Marking after removal of tumor
- Marking of suspicious lesions
- Tumor marking during neoadjuvant chemotherapy
- For marking the biopsy site
- Orientation for planning radiation therapy

Advantages at a glance

- Shape and size adapted for the application within the MRI gantry
- 3D marker design ensures firm anchoring in the tissue
- Optimal visibility in all positions
- Cannula is MRI compatible
- Marker is MRI compatible up to 3 Tesla
- Marker is approved as implantable material (Nitinol)
- Ergonomic handle for single-handed operation
- The marker is supplied preloaded

Product	Order No.	Cannula size
Tumark® MRI	601 570	18G/1,2 x 120 mm


MRI Tuloc®

Stability

The MRI Tuloc® Localization System serves for the pre-operative marking of suspicious tissue under MRI control. Thanks to increased pressure stability and an extremely sharp bevelling of the cannula and the wire tips, MRI Tuloc® even allows for an easy and safe penetration of solid tumor tissue. If a repositioning is required, the wire can be easily pulled back into the cannula for a correction of the position. The wire of the Tuloc® has markings for secure location for the user. Two proximal markings indicate whether the arches of the wire are inside or outside of the cannula. The distal markings at distances of 2 cm give information on the distance to the tumor.



Beveled arches facilitate penetration



Laser-cut double arch

Monofilament wire provides for best possible form and pressure stability

This marking indicates that the arches are still inside the cannula

Markings every 2 cm provide information on the distance to the distal tip of the wire

This marking shows that the arches are completely unfolded

4 cm to the distal tip of the wire

2 cm to the distal tip of the wire

Advantages at a glance


- Wire and cannula are MRI compatible
- Repositionable localization wire
- High rigidity and stability of the cannula
- Monofilament wire provides for best possible form and pressure stability
- Distal depth markings show information on the distance to the tumor
- Extreme sharpness for precise and atraumatic puncture
- Optimal visibility in the MRI and good palpability during the operation

MRI Tuloc®			
Order No.	Gauge	Diameter	Length
601 649	20 G	0,95 mm	90 mm
601 651	20 G	0,95 mm	120 mm

MRI Duo-System


Flexible

The SOMATEX® MRI DUO System is a correctable localization system for the pre-operative marking of non-palpable, suspicious mammary lesions under MRI control. If a repositioning is required, the localization wire can be easily pulled back into the cannula and reused for a correction of the position.



Fixation piece

Repositionable double arches



Advantages at a glance

- Wire and cannula are MRI compatible
- High rigidity and stability of the cannula
- Correctable and flexible marking wire
- Sharp wire tips for firm lesions


Varying threading of the wires every 50mm allows for depth orientation



Repositionable and flexible localization wire

Marking to identify MRI version

Sharp needle tip



Puncture needle with high stability

MRI Duo-System			
Order No.	Gauge	Diameter	Length
601 607	20 G	0,95 mm	50 mm
601 609	20 G	0,95 mm	90 mm
601 611	20 G	0,95 mm	120 mm

MRI Biopsy Handy

Lightweight

SOMATEX® **MRI Biopsy Handy** is a semi-automatic device for obtaining histologically valuable tissue material from a variety of soft tissues and organs. The novel material and special alloy allow for optimal visibility in the MRI and therewith the exact positioning of the specimen notch. The tissue material is separated by means of an extremely fast forward movement of the exterior cutting cannula. When complemented with the MRI coaxial puncture needle, one Biopsy Handy can be reused for multiple tissue sampling.

MRI Biopsy Handy			
Order No.	Gauge	Diameter	Length
603 108	18 G	1,20 mm	100 mm
603 110	18 G	1,20 mm	150 mm
603 112	18 G	1,20 mm	200 mm
603 118	16 G	1,60 mm	100 mm
603 120	16 G	1,60 mm	150 mm
603 122	16 G	1,60 mm	200 mm
603 128	14 G	2,00 mm	100 mm
603 130	14 G	2,00 mm	150 mm
603 132	14 G	2,00 mm	200 mm

Matching MRI coaxial puncture needle

MRI Coaxial Puncture Needle			
Order No.	Gauge	Diameter	Length
601 400	18 G	1,20 mm	100 mm
601 405	18 G	1,20 mm	150 mm
601 408	16 G	1,60 mm	45 mm
601 410	16 G	1,60 mm	95 mm
601 412	16 G	1,60 mm	145 mm
601 418	15 G	2,00 mm	43 mm
601 420	15 G	2,00 mm	93 mm
601 422	15 G	2,00 mm	143 mm
601 428	13 G	2,40 mm	41 mm
601 430	13 G	2,40 mm	91 mm
601 432	13 G	2,40 mm	141 mm

Matching MRI Biopsy-Handy

MRI Coaxial

Puncture Needle

One for all

The **MRI coaxial puncture needles** are multipurpose cannulas and can be used for initial punctures, as sheaths to guide wires, and as puncture sheaths for both the SOMATEX® MRI Biopsy Handy and the standard version Biopsy Handy. The novel material, special alloy, and particularly sharp trocar tip allow for optimal visibility in the magnetic resonance imaging, at maximum stability, and minimum puncture trauma.



Application

- Initial puncture
- Coaxial cannula for biopsy instruments
- Guide needle for catheter and guide wires
- Injection and drainage cannula

MRI Chiba Needle

Universal

The SOMATEX® **MRI Chiba Needle** is a multipurpose cannula. It is used for the injection of drugs and for fine needle aspiration biopsies. The novel material, special alloy, and the particularly sharp ultra cut allow for optimal visibility in the MRI, at maximum stability, and minimal puncture trauma.


Application

- Pain therapy
- Fine needle aspiration biopsy

Advantages at a glance



- MRI compatible
- High rigidity and stability
- Minimal puncture trauma due to particularly Sharp ultra cuts
- Depth stopper

MRI Chiba Needle			
Order No.	Gauge	Diameter	Length
601 250	22 G	0,70 mm	100 mm
601 255	22 G	0,70 mm	150 mm
601 260	20 G	0,95 mm	100 mm
601 265	20 G	0,95 mm	150 mm



Advantages at a glance

- MRI compatible
- High strength and stability
- Optimal visibility in the MRI
- Easy handling with one hand
- Light weight
- Sharp cut of stylet and cannula tip
- Optionally available with coaxial sheath



SOMATEX[®]

MEDICAL TECHNOLOGIES GMBH

Rheinstr. 7d
D-14513 Teltow
Germany

Tel. +49 (0) 3328 30 76-0
Fax +49 (0) 3328 30 76-99
info@somatex.com
www.somatex.com