

### \* Biopsy Device Compatibility<sup>#</sup>

| Biopsy Devices                                      | MAM3001 | MAM3002 | MRM4002 | MRM4008 |
|-----------------------------------------------------|---------|---------|---------|---------|
| ATEC 0909-20, 9-Gauge Standard                      | ✓       | ✓       | ✓       |         |
| ATEC 0909-12, 9-Gauge Petite                        | ✓       | ✓       |         |         |
| ATEC 0912-20, 9-Gauge Long                          | ✓       | ✓       | ✓       |         |
| ATEC 0912-12, 9-Gauge Long Petite                   | ✓       | ✓       |         |         |
| ILS 0914-20, 9-Gauge Introducer Localization Set    |         |         | ✓       |         |
| ILS 0914-12, 9-Gauge Introducer Localization Set    |         |         | ✓       |         |
| EnCor® 7-Gauge CP01-7G                              |         |         |         | ✓       |
| EnCor® 10-Gauge CP01-10G                            |         |         | ✓       |         |
| EnCor® MRI Probe & Introducer Set ECPMRPAX-10-Gauge |         |         | ✓       |         |
| Coaxial Cannula for 10-Gauge x 118mm (VC10118)      | ✓       | ✓       |         |         |
| Coaxial Cannula for 10-Gauge x 140mm (VC10140)      |         |         | ✓       |         |

| Code    | Description (clip)                                                      | How supplied |                                                                                       |
|---------|-------------------------------------------------------------------------|--------------|---------------------------------------------------------------------------------------|
| MAM3001 | MammoMARK 11-Gauge                                                      | 10 per box   |  |
| MAM3002 | MammoMARK2 11-Gauge<br>(differentiated titanium marker for second site) | 10 per box   |  |
| MAM3008 | MammoMARK 8-Gauge                                                       | 10 per box   |  |
| MRM4002 | MammoMARK2 11-Gauge MR                                                  | 10 per box   |  |
| MRM4008 | MammoMARK2 8-Gauge MR                                                   | 10 per box   |  |
| MAM3014 | CorMARK 14-Gauge                                                        | 10 per box   |  |

For complete product details see Instructions for Use.

CorMARK and MammoMARK are trademarks of Devicor Medical Products, Inc..

<sup>#</sup>Data on file at Devicor Medical Products, Inc.