

Biopsy Valves

Single Use

Biopsy Valves

- Covers and Seals the Opening of the Accessory Channel Inlet
Supports the Insufflation and Suction Function
Avoids the Leakage of Potentially Hazardous Bodily Fluids
- Versions with Irrigation Line and Luer-Lock Adapter
Connectivity to a Syringe or External Water Pump



Biopsy Valves

Single Use



Biopsy Valves

The latex-free biopsy valves are used to cover and seal the opening of accessory channel inlet of gastroscopes, colonoscopes and duodenoscopes. The valves:



- a) facilitate the insufflation and suction function during an endoscopic procedure.
- b) avoid the leakage of potentially hazardous bodily fluids for the protection of the patient and medical staff, for example while changing endoscopic instruments.

The version, featuring an irrigation line allows to connect a syringe or an external water pump.

Biopsy Valves – single use, individually packed in a dispenser

Catalog No.	Description	Colour	Unit	Product
GAC-04-00-001	Biopsy Valve, compatible to endoscopes of: → Olympus → Fujinon (500 Series and newer) → Storz	Blue	Box of 100	
GAC-04-01-001	Biopsy Valve, compatible to endoscopes of: → Pentax	Red	Box of 100	

Biopsy Valves with Irrigation Line and Luer-Lock Adapter – single use, individually packed in a dispenser

Catalog No.	Description	Colour	Unit	Product
GAC-04-00-002	Biopsy Valve, compatible to endoscopes of: → Olympus → Fujinon (500 Series and newer) → Storz	Blue	Box of 50	
GAC-04-01-002	Biopsy Valve, compatible to endoscopes of: → Pentax	Red	Box of 50	

Ver. 1.0 (2011-01)