



sensitive. versatile. confident.  
**just like you.**



**Neoprobe**  
Gamma Detection System

# Proceed with confidence.



## Easy System Set-Up

- + No need for calibration
- + Integrated Bluetooth® wireless technology

## Options

- + Dynamic Pitch – lymphatic mapping
- + Binary Pitch – parathyroid and FDG isotope procedures
- + Autoranging – Focus audible feedback on detected energy level for more precise localization

## Software Upgradeable

- + Provides peak performance in changing surgical procedures

## Bluetooth Wireless Probes

- + Cordless, ergonomic probes, angled and straight
- + Laparoscopic cordless probe
- + Internally collimated 9mm straight probe
- + More sensitive than the NPR14 corded probe
- + High Energy (F-18) Probe detects radioactive isotope Fluorine-18, called 18F-FDG, and other high-energy radionuclides

## User Selected Threshold

- + Sets windows for most common radioisotopes used in surgical applications:  $^{125}\text{I}$ ,  $^{57}\text{Co}$ ,  $^{99\text{m}}\text{Tc}$ ,  $^{111}\text{In}$ ,  $^{131}\text{I}$ ,  $^{18}\text{F}$
- + A Seventh OPEN Window. Dual isotope mode for simultaneous scanning of  $^{99\text{m}}\text{Tc}$  and  $^{125}\text{I}$

## Ordering Information

Neoprobe™ GDS Control Unit with Software Included (Requires power cord, sold separately)	NPCU3
9mm Reusable Probe with Bluetooth® II Technology	NPB09S
11mm Reusable Bluetooth® Laparoscopic Probe (No additional cable required. For use with Bluetooth® cable system only)	NPB11L
14mm Reusable Bluetooth® Probe, Angled	NPB14A
14mm Reusable Bluetooth® Probe, Straight	NPB14S
Neoprobe™ High Energy Reusable Probe (Requires High Energy Probe Cable (NPAF18), sold separately. For use with Neoprobe™ Control Unit, model 2200 and 2300 only)	NPRF18

For more information about Mammotome®, Neoprobe™ and other products, or to place an order, contact your Mammotome Sales Representative or call 1-877-9-A-MAMMO (1-877-926-2666).



Neoprobe™ and neo2000™ are trademarks of Devicor Medical Products, Inc. Bluetooth® word mark and logos are owned by Bluetooth® SIG, Inc. and any use by Devicor Medical Products, Inc. or its representatives, is under license. For complete product details, see Instructions for Use. ©2012 Devicor Medical Products, Inc. All rights reserved.